


## NeuroNexus Adaptors

The picture below illustrates a Plexon NeuroNexus adaptor, which connects a Plexon headstage to a NeuroNexus Probe.




Plexon offers a variety of adaptors to interface our headstages to acute NeuroNexus Probes. Please contact our [sales department](#) to discuss the appropriate configuration for your probe.

## NeuroNexus Probes


### Adaptors for acute NeuroNexus Probes

**Description:** This is a top and bottom view of an adaptor to interface two 8-channel Plexon headstages to a 16-channel NeuroNexus Probe. (Other configurations available)

Top



Bottom



### Remapper

**Description:** Plexon also offers a solution to map the channels of the NeuroNexus Probes to consecutive Plexon preamp channels. Contact our [sales department](#) to determine which remapper is compatible with your probe.



## Adaptors for Chronic 16 and 32-channel NeuroNexus Probes

No adaptors are required to interface to the NeuroNexus chronic probes. The Plexon 16V and 32V headstages interface directly to the Omnetics connector on the probe.